

Ceramic Spittoon



Four Way Foot Control



Modular Control



- ❖ Semi Electrical Dental Chair with Chair Mounted Dental Unit.
- ❖ Electrically controlled Up/Down movements.
- ❖ Manually operated Backrest.
- ❖ Shadow less operating light with dual intensities (9V/12V).
- ❖ Ceramic spittoon.
- ❖ Air ventury HV Suction unit.
- ❖ Horizontally adjustable instrument tray with 3 point modular control for 3-way syringe (with additional tip) & Airtor control Plus Micro motor complete (without Hand pieces).
- ❖ Pneumatic Dental operator's stool.



Colours and features shown may vary.