

# THE *Ergodental* DESIGNER FURNITURE


Novo Sar

## YOUNG 225

### WORK BENCH

Work benches are arranged for the fitting up of the wide range of accessories available.

Double workbench



## YOUNG 203

Single workbench



### Light Kit



ISY 210 - Light Kit with 36W  
neon tube & shelf for Young 203

ISY 230 - Light Kit with 58W  
neon tube & shelf for Young 225

### THE STANDARD VERSION INCLUDES

- Retractable air gun
- Pressure regulation and manometer
- Fuse and housing for micromotor with power socket.
- Gas tap

### OPTIONALS

Dust extraction with vacuum pump and filter drawer

Wide range of extraction peg

Drawer front with several aesthetics

and functional solutions

### Dust Extraction Peg



Colours and features shown may vary.

### Corporate Office:

428, Laxmi Plaza, Off New Link Rd., Andheri (W), Mumbai- 400 053. (India)  
Tel : +91 - 22 - 4017 5555 • Fax: 6697 6191  
E-mail : info@novodental.co.in , www.novodental.co.in



**NOVODENT**  
EQUIPMENTS & MATERIALS PVT. LTD.

### Production Unit :

C- 33 / 35, Industrial Estate, Thattanchavady, Pondicherry- 605 009.  
Tel : +91 - 413 - 420 2171 / 72 / 91 / 420 2271 • Fax: +91 - 413 - 2245614  
E-mail : nempl@novodental.co.in / nempl@gmail.com