

TRULY SPARKLED

Glass Spittoon



Six Way
Foot Control



Modular
Control



Swiveling
Assistant Console



Available
Upholstery Colour Samples

Colours and features shown may vary.

- ◆ Fully electrical dental chair
- ◆ Feathertouch control with four working programs
- ◆ Shadowless operating light with dual intensities (9V/12V)

- ◆ Coloured Glass Spittoon
- ◆ Motorized HV/LV suction with auto drain
- ◆ Fully adjustable instrument tray with 4 point modular control having 3-way syringe, two coupling for Airtor, Micro motor complete (without any Hand pieces)

- ◆ Stainless steel instrument tray and LED X-ray viewer
- ◆ Six way foot control
- ◆ Swiveling assistant 3-point console
- ◆ Pneumatic Dental operator's stool

Scintila Continental

TRULY SPARKLED

Glass Spittoon



Six Way
Foot Control



Modular
Control



Swiveling
Assistant Console



Swiveling
Assistant Console



Available
Upholstery Colour Samples

Colours and features shown may vary.

- ◆ Fully electrical dental chair
- ◆ Feathertouch control with four working programs
- ◆ Shadowless operating light with dual intensities (9V/12V)

- ◆ Coloured Glass Spittoon
- ◆ Motorized HV/LV suction with auto drain
- ◆ Fully adjustable instrument tray with 4 point modular control having over head delivery system, 3-way syringe, two coupling for Airtor, Micro motor complete (without ANY Hand pieces).

- ◆ Six way foot control
- ◆ Swiveling assistant 3-point console
- ◆ Pneumatic Dental operator's stool